PHARMA, FOOD & BIOTECH PRODUCTS







simultaneous CNC machines with CAD connection, were necessary to achieve such a performance capability. These machines produce GEMÜ M600 multi-port valves where all the valve weirs and cross connections are integrated "into one component".

These features save our customers costs and enable us to guarantee the highest possible service reliability.



GEMÜ - "THE ENGINEERS CHOICE"

Diaphragm and globe valves in stainless steel
Autoclavable, aseptic, almost deadleg-free
Surface finishes down to Ra 0.25 µm, electropolished
Straight through and multiport valves, valve manifolds
EPDM, FPM, PTFE diaphragms and seals

Manual, pneumatic or motorized actuators

Nominal sizes DN 2 - 100

Measurement, control and regulating devices for pressure, temperature, level and flow rate Integrated field bus connections













